

HIV SEROCONVERSION PANELS

Zeptomatrix has the most comprehensive range of Seroconversion panels in the blood screening industry. This is one of the reasons that Zeptomatrix panels are the first choice of many International Regulatory Agencies and Reference Laboratories.



- **Sterile, unadulterated, therapeutic grade source plasma from asymptomatic untreated donors (frozen within 30 minutes of collection)**
- **Maximum 2 freeze thaw cycles on any panel member**
- **Complete full Seroconversion profiles from negative to antibody positive**

NEW PANELS!

PANEL NUMBER	DONOR NUMBER
HIV9096	1043487

Anti-HIV-1 False Positive

PANEL NUMBER	DONOR NUMBER	PATIENT DRAWS	DAYS SPANNED	HIV RNA	HIV-1 bDNA	HIV Ag	ANTI-HIV	3RD GEN	4TH GEN
								HIV AB 1/2+0	HIV-Ag/Ab COMBO
HIV9027	65633	9	29	*	POS	NEG	POS	POS	POS
HIV9078	79065	5	70	*	*	POS	NEG	NEG	NEG

Anti-HIV-1

PANEL NUMBER	DONOR NUMBER	PATIENT DRAWS	DAYS SPANNED	HIV RNA	HIV-1 bDNA	HIV Ag	ANTI-HIV	HIV AB 1/2+0	HIV-Ag/Ab COMBO
HIV6243	62238	10	35	POS	*	POS	POS	POS	POS
HIV6244	61694-74	14	34	POS	*	POS	POS	POS	POS
HIV6247	63602	10	31	POS	POS	POS	POS	POS	POS
HIV6248	63331	7	26	POS	POS	POS	POS	POS	POS

**Do you have any question or would you be interested to receive our quotation?
Please contact us by e-mail (sales@helveticahealthcare.com)
or by phone (+41 22 534 94 54)**

PANEL NUMBER	DONOR NUMBER	PATIENT DRAWS	DAYS SPANNED	HIV RNA	HIV-1 bDNA	HIV Ag	ANTI-HIV	HIV AB 1/2+0	HIV-Ag/Ab COMBO
HIV9010	65376	8	50	*	POS	POS	POS	*	*
HIV9011	64954	11	39	*	POS	POS	POS	POS	POS
HIV9012	65389	8	23	*	POS	POS	POS	POS	POS
HIV9013	65404	7	26	*	POS	POS	NEG	NEG	POS
HIV9014	65522	7	31	*	POS	POS	POS	POS	POS
HIV9015	65661	8	38	*	POS	POS	POS	POS	POS
HIV9016	65790	10	35	*	POS	POS	NEG	NEG	POS
HIV9017	65907	6	36	*	POS	POS	POS	POS	POS
HIV9018	66575	11	36	*	POS	POS	POS	POS	POS
HIV9019	65685	3	38	*	POS	POS	POS	POS	POS
HIV9020	67512	22	98	*	POS	POS	POS	POS	POS
HIV9021	67485	17	58	*	POS	POS	POS	POS	POS
HIV9022	64578	9	33	*	POS	POS	POS	POS	POS
HIV9023	67706	22	86	*	POS	POS	POS	NEG	POS
HIV9024	64469	12	54	*	POS	POS	POS	NEG	POS
HIV9025	67996	12	92	*	POS	POS	POS	POS	POS
HIV9026	68205	7	45	*	POS	POS	POS	POS	POS
HIV9028	67860	7	56	*	POS	POS	NEG	NEG	POS
HIV9029	68262	12	46	*	POS	POS	POS	NEG	POS
HIV9030	68582	16	55	*	POS	POS	POS	POS	POS
HIV9031	66048	19	158	*	POS	POS	POS	POS	POS
HIV9032	68106	14	56	*	POS	POS	POS	POS	POS
HIV9033	66686	16	85	*	POS	POS	POS	POS	POS
HIV9034	66632	13	54	*	POS	POS	POS	POS	POS
HIV9075	62216	5	48	POS	*	POS	POS	POS	POS
HIV9076	63753	9	74	POS	*	POS	POS	POS	POS
HIV9077	75018	27	92	*	*	POS	POS	POS	POS
HIV9079	75062	25	98	POS	*	POS	POS	POS	POS
HIV9081	63215	4	34	POS	*	POS	POS	POS	POS
HIV9082	77840	4	21	POS	*	POS	POS	*	*
HIV9084	77600	4	49	POS	POS	NEG	POS	*	*
HIV9089	65376	6	24	*	POS	POS	POS	POS	POS
HIV9096	1043487	6	22	POS	*	*	POS	*	*
HIV12007	73695	9	134	POS	*	NEG	POS	POS	POS
HIV12008	73698	13	43	POS	*	POS	POS	POS	POS

All products are for Research Use Only and not for in vitro Diagnostic Use

Package Inserts are provided for product evaluation purposes only and are not intended to be used in place of the Package Insert shipped with the product.

[Link to the webpage](#)

**Do you have any question or would you be interested to receive our quotation?
Please contact us by e-mail (sales@helveticahealthcare.com)
or by phone (+41 22 534 94 54)**